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## GB Instructions for use

### Elastomeric bite registration material on vinyl polysiloxane basis, addition-curing.

Thank you for the confidence you have shown in our company. You have made the right decision with the purchase of this bite registration material. Please be sure to read these instructions for use through carefully before starting to use the product! The bite registration materials are only intended for use in dental applications by correspondingly qualified professionals. The medical device must be employed responsibly and with due caution.

#### Product description:

All bite registration materials are produced from high-quality raw materials and suitable for a wide range of applications thanks to their different final hardnesses. The materials offer very high thixotropism, are easy to apply thanks to the automatic safety cartridge system, can be removed from the mouth with ease and display long-term dimensional stability. The USP of DC Bite extra hard scan is the optical reflection activity, which makes it especially suitable for powder-free optical 3D registration and data capture in the field of CAD/CAM/CIM.

#### Indication:

DC Bite NF + DC Bite Apple: Production of long-lasting, reusable bite registrations; fixation of support and facebow registrations. DC Bite extra hart scan: Optical data capture in the field of CAD / CAM / CIM (antagonist registration); standard bite registration.

#### Bite registration material in 50 ml safety cartridges

Insert the cartridge in the 1:1 mixing gun, remove the end cap and dispense material until material is flowing evenly from both openings. Now place the mixing tip (single-use) over the openings (ensure that the colour code of the end cap and mixing tip are identical) and ensure that it is positioned correctly. To close it, turn the tip 90° clockwise. If necessary, use an intraoral tip (single-use). Operate the mixing gun again and dispense material in the desired amount. The used mixing tip remains on the cartridge and serves as a cap. Prior to the next use, remove the used mixing tip, check the opening for a plug that forms in

#### Technical data and product features

Product	Colour	Dosage (1:1)	Working time incl. mixing time*	Time in mouth	Total setting time	Hardness (+/-5)	Linear dimensional change
DC Bite NF	orange	automatic	30 sec.	30 sec.	1 min.	Shore A90	0,1 %
DC Bite Apple	green	automatic	30 sec.	60 sec.	90 sec.	Shore A95	0,1 %
DC Bite extra hart scan	yellow	automatic	30 sec.	30 sec.	1 min.	Shore D45	0,1 %

\* At 23°C / 73.4°F. The times specified in the table above are based on 23°C and a relative humidity of 50 +/- 10%. As a general rule: High temperatures accelerate setting and lower temperatures delay it.

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#### Storage:

Store in a dry place at between 5°C and 27°C (41°F and 80.6°F).

#### Shelf life:

See date printed on package / label. Do not use after expiry date.

#### Processing temperature:

Room temperature between 18°C and 25°C (64.4°F to 77°F).

#### Contraindication:

The product must not be used in patients with known allergies to flavourings or other ingredients (see safety data sheet or ask the manufacturer). This should be clarified before use!

#### Warnings:

##### Eye contact:

Avoid contact with eyes. In case of contact, rinse eyes immediately with plenty of clean water and consult an ophthalmologist without delay.

##### Swallowing:

If swallowed, rinse out the mouth and drink plenty of water. Always contact a doctor.

#### Allergic reactions:

It is not possible to rule out allergic reactions in hypersensitive people. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately and a doctor should be consulted. Should a serious incident occur unexpectedly in connection with this medical device, it must be reported to the manufacturer and the competent authority in your Member State.

#### Please also observe the warnings in the safety data sheet!

Information at [www.dental-central.de](http://www.dental-central.de)  
Errors and omissions excepted

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