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Patientenzielgruppe

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en Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended users

Dentist

Composition

Qualitative composition: 1,4-butandiol-dimethacrylate, Bis-GMA, aliphatic urethane dimethacrylate, butylated hydroxytoluene, camphorquinone, silicon dioxide, glass ceramic, iron oxide pigments, 2-(2H-Benzotriazol-2-yl)-p-cresol, 2-ethylhexyl-4-(di-methyl-amino) benzoate, diphenyl (2,4,6-trimethyl-benzoyl) phosphine oxide

Quantitative composition: 1,4-butandiol dimethacrylate 5%, aliphatic urethane dimethacrylate 22 %, Bis-GMA 18 %, silicon dioxide 6%, activators and UV absorber <1 %, glass ceramic 48%, pigments <1 %

Intended purpose of the product

DC FlowFill is a light-curing, radiopaque, free-flowing composite filling material.

Indications

Cavities of classes III, IV and V.

Contraindications

Do not use in patients with known allergy to the ingredients.

Performance characteristics of the product

DC FlowFill is a free-flowing composite that cures with light.

Specific precautions

None

Application

1. Clean and dry the preparation. Determine the colour.
2. Prepare the cavity according to the general rules of adhesive technique. Bevel the edges of the enamel.
3. Apply etching gel to the enamel edges for 20 to 35 seconds.
4. Completely rinse off the etching gel.
5. Dry with oil-free air. Avoid over-drying the dentin.
6. Apply adhesive as described in the usage instructions.
7. Select the syringe in the desired colour. Remove the syringe cap. Screw on the Luer Lock canula.
8. DC FlowFill is applied directly by pressing the syringe. Apply DC FlowFill in layers no thicker than 2 mm to ensure complete curing.
9. Leave each layer to light cure for 30 to 40 seconds.
10. Repeat steps 8 and 9 until the cavity is filled.
11. Remove excess, contour and finish with diamond finishers and flexible discs of decreasing grain size.

Warnings

Residual risks

Please observe carefully

Light units should be emitting at 450 nm at 500 mW/cm2 and should be checked regularly. Place the light as close as possible to the filling material. Do not remove the oxygen-inhibited layer, as this is needed for bonding with the next layer. Using a die is recommended. Follow the instructions for use of other materials used with DC FlowFill.

Safety information

H317 May cause allergic skin reactions

H318 Causes serious eye damage

Use is contraindicated if drainage or the prescribed application method is not possible. Do not use directly on opened pulp or if there is a known allergy to any ingredient. No side effects have been reported to date. In individual cases, hypersensitivity to components of the material may occur. Avoid skin contact with the paste. In case of accidental contact, wash immediately with soap and water. In case of eye contact, rinse immediately with plenty of water and consult a doctor. Avoid materials containing eugenol, moisture and oily air, as they can prevent polymerisation at the point of contact. Discolouration may occur when using cationic mouthwashes, plaque disclosers or chlorhexidine.

Side effects

Adverse side effects are not known. Risks in case of overdose are not known.

Interactions

Interaction of the product and its metabolites with other medical devices, medicinal products or other substances is not known.

Storage

Dry, do not store above 22 °C.

Minimum shelf life

Do not use after the expiry date.

Variants

[REF 716321](#) 1x syringe A1 with 1,0ml + cannulas

[REF 716322](#) 1x syringe A2 with 1,0ml + cannulas

[REF 716323](#) 1x syringe A3 with 1,0ml + cannulas

[REF 716324](#) 1x syringe A3,5 with 1,0ml + cannulas

[REF 716327](#) 1x syringe B2 with 1,0ml + cannulas

Disposal

Unused material can be disposed of with household waste after curing.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

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