# DC Glas-IonomerCem

DC DentalCentral GmbH • Owidenfeldstraße 6 • D-30559 Hannover • Germany • www.dental-central.de

Stand 02 / 2021

## D Gebrauchsinformation

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## GB Instruction for use

DC Glas-IonomerCem is an X-ray opaque glassionomer cement that has been specifically developed as an underfilling material underneath composite and amalgam fillings. DC Glas-IonomerCem increases the biocompatibility of composites and improves the border seal. Adhesion to the dentin and enamel is by way of the reaction of the carboxylate groups of polyacrylic acid with the calcium or hydroxylapatite and collagen. This direct chemical bonding ensures that far more healthy dental material is saved because undercutting is not necessary.

**Classification**  
DIN EN ISO 9917-1:2004-04 requirements are fulfilled.

**Preparing the cavities**  
A maximum of healthy dentine must be retained when preparing the cavities. The cavity must be shaped in such a manner that the cement layer is not thinner than 0.5 mm at any point.

**Pulp protection**  
Rinse the prepared cavities thoroughly with water and dry with oil-free blown air. Use calcium hydroxide preparations as a pulp-capping for very deep cavities.

**Preparing the dentin**  
Clean the freshly prepared dentin surface, rinse thoroughly with water and dry with blown air. Avoid dehydration.

**Mixing DC Glas-IonomerCem**  
Use one level dosing scoop for two full drops of liquid (without air bubbles); hold the bottle vertically downwards and apply light pressure to squeeze out one drop. First mix only half the amount of powder, then continue with the other half of the remaining powder until the required consistency is achieved. Use only a small area of the mixing block in order to minimize water evaporation.

**Application of DC Glas-IonomerCem**  
The mixed material is lightly pressed into the cavity. DC Glas-IonomerCem is workable for 1'40 minutes (incl. mixing). The material should not be handled any more when hardening sets in. The hardening time is 3 - 3'40 minutes. Only cover the cavity base with DC Glas-IonomerCem if amalgam is to be used. The amalgam can be applied after 3'40 minutes. The dentin should be completely covered with DC Glas-IonomerCem when composites are used as filling material. Treat the entire cavity surface (incl. the cement

surface) with a dentin/enamel adhesive (e.g. DC Dentine-Composite Bond) according to the manufacturer's instructions. The composite filling is applied according to the manufacturer's instructions.

**Side-effects**  
Side-effects with this medical product are extremely rare if it is properly prepared and applied. Should you encounter undesirable side-effects - even in case of doubt - then please contact us and give a precise description of the circumstances and symptoms.

**Processing instructions**  
Ensure that the processing instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the circumstances.

**Storage**  
Shelf life: See label! Opened packs can be kept for 6 months provided they are kept dry and stored at room temperature.

**WARNING!** Do not allow non-bonded material to come into contact with the eyes. Should this nevertheless inadvertently happen, then rinse immediately with plenty of water.

**Please note**  
Instruments must be immediately cleaned under flowing water after they have been used.

**Packaging**  
15 g powder / 10 g liquid + dosing scoop

*The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.*

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